

ENVIRO-CLEAN® Universal Cartridges

Compatible with the Horizon SPE-DEX® 4790 Automated Extraction System



SPE-DEX® is a registered trademark of Horizon Technology, Inc.

The Enviro-Clean® Universal Cartridge is the choice of modern contract labs because it is inexpensive and easy to use. Enviro-Clean® sorbents, UCT polypropylene, and PTFE frits reproducibly offer clean blanks, consistent extractions, and consistent flow rates. Designed for the environmental lab, the cartridge is made to handle large volumes of waste water. An optional bottle holder is available for continuous feed from Boston Round and wide mouth bottles.



amchro GmbH • Tel: 06190-978 117 • Fax: 06190-978 107 • www.amchro.com

ENVIRO-CLEAN® Universal Cartridges

The cartridge will fit all standard manifolds and disk manifolds with adapter.

Product Name	Part Number	Sorbent Amount/ Tube Volume	Units	Description / Application
UNIVERSAL C18	ECUNIC18	1100mg/83mL	8	1100 mg of endcapped C18 for pesticides, PCBs and a large assortment of applications.
UNIVERSAL 525	ECUNI525	1500mg/83mL	8	1500 mg of our special C18 blend. This cartridge is specifically designed for EPA Method 525.
UNIVERSAL PAH	ECUNIPAH	2000mg/83mL	8	2000 mg of C18 specifically designed for PAH extraction.
UNIVERSAL DRO	ECUNIDRO	1500mg/83mL	8	1500 mg of C18 specifically designed for the extraction of diesel range organics and similar compounds.
UNIVERSAL OIL & GREASE	ECUNIOAG	4000mg/83mL	15	4000 mg of large particle C18 with an assortment of PE frit filters. No more liquid/liquid emulsions or clogged disks.
UNIVERSAL 521	ECUNI521	2000mg/83mL	8	2000 mg of our special blend of activated carbon; Specifically designed for EPA Method 521 "Nitrosamines in Drinking Water" and other methods where activated carbon is the sorbent of choice.

