


A 3-IN-1 SOLUTION



Universal Vacuum Manifold




UNITED CHEMICAL TECHNOLOGIES, INC.

2731 Bartram Road,
Bristol, PA 19007

PHONE: (800) 385-3153 • FAX : (215) 785-1226 • WEB: www.unitedchem.com

United Chemical Technologies, Inc.
Universal Vacuum Manifold

Accommodates the following configurations:




24 Well Plate

- Extraction plate contains (24) 6 mL SPE columns
- UCT SPE columns can be packed in 500 mg, 1 and 2 gram formats.
- Collect samples directly into microplate or analytical vials.



48 Well Plate

- Extraction plate contains (48) 4.5 mL SPE columns in a staggered array
- Elute samples directly into 12 x 75 mm test tube or 2 mL analytical vials.



96 Deep Well Plate

- Accommodates standard 96 deep well plate - 2 mL wells
- Collect directly into microplate



2731 Bartram Road,
Bristol, PA 19007